

Merivaara Emergo Patient Trolley



The multipurpose Emergo patient trolleys (stretchers) are used for patient transportation, emergency, recovery, ICU and minor procedures. The safe and sturdy construction makes it ideal for heavy use.

The versatile range of options and adjustments ensure efficient patient transfer, patient comfort and ergonomic working positions for health-care staff. The adjustable mattress base provides comfortable positions for patients to recover eliminating unnecessary patient transfer from the trolley (e.g. in sitting position).

Large, easy rolling castors provide excellent mobility and the optional fifth wheel improves overall mobility reducing the strain for the staff. Emergo has many attachment points for the IV-rod and other accessories. The wide range of accessories available will serve the demanding needs of nursing staff for many years.

Main features:

1. Modular structure that allows several different configurations
2. Clear and easily controllable functions
3. It is easy and light to carry out different care positions

The range:

- 6250: 2-part bed section with hydraulic height adjustment and gas spring assisted Trendelenburg and back section
- 6270: 2-part X-ray bed section with hydraulic height adjustment and gas spring assisted Trendelenburg and back section
- 6280: 2-part X-ray bed section with subframe, hydraulic height adjustment and gas spring assisted Trendelenburg and back section
- 6350: 3-part bed section with hydraulic height adjustment and gas spring assisted Trendelenburg, back section and leg section
- 6370: 3-part X-ray bed section with hydraulic height adjustment and gas spring assisted Trendelenburg, back section and leg section

Merivaara Emergo Patient Trolley

Model	6250	6270	6280	6350	6370
Base sections	2	2	2	3	3
X-ray platform	No	Yes	Yes	No	Yes
Subframe for equipment	No	No	Yes	No	No
Width	750 mm	750 mm	750 mm	750 mm	750 mm
Patient surface width	705 mm	705 mm	705 mm	705 mm	705 mm
Length	2075 mm	2075 mm	2150 mm	2075 mm	2075 mm
Recommended mattress size	650 x 1950 mm	650 x 1950 mm	650 x 1950 mm	650 x 1950 mm	650 x 1950 mm
Height adjustment	560 - 930 mm	625 - 995 mm	625 - 995 mm	560 - 930 mm	625 - 995 mm
Safe Working Load (SWL)	225 kg	225 kg	225 kg	225 kg	225 kg
Trendelenburg	14°	14°	14°	14°	14°
Anti-Trendelenburg	7°	7°	7°	7°	7°
Back section adjustment	0 - 70°	0 - 70°	0 - 70°	0 - 70°	0 - 70°
Leg section adjustment	-	-	-	0 - 45°	0 - 45°
Castors	Central locking, antistatic 150 mm	Central locking, antistatic 150 mm	Central locking, antistatic 150 mm	Central locking, antistatic 150 mm	Central locking, antistatic 150 mm
Country of manufacture	Finland	Finland	Finland	Finland	Finland
CE marked	Yes	Yes	Yes	Yes	Yes

Other features: Yellow bumpers in each corner | pedals for height adjustment on both sides

Factory installed options:

Wire basket 582 under back section | wire basket 583 under foot section | xygen bottle holder under the mattress base 2 or 5 kg | paper roll holder | central locking castors 200 mm XSP and 5th wheel | ABS storage boxes small or large | fixed foldable infusion rod for subframe

Other accessories:

Arm rest | infusion rod 4 hooks | accessory rail, pair | accessory basket | restraining straps | infusion rod holder | trolley end | side rails, chromed with yellow padding | foldable push bar, chromed with yellow padding | D-pillow, visco-elastic, including fixing straps 21x49x10cm

Mattresses:

Mattress 50 mm, antistatic, black, side fixing, bottom velcros

Mattress 80 mm, antistatic, black, side fixing, lifting handles, bottom velcros

Mattress 100 mm VEF antistatic, black, side fixing, lifting handles, bottom velcros